


APLW 051

Injector nozzle

with plastic support, star-shaped with clamp, Ø 4, length 110 mm.

EAN: 4002516116127 / Material number: 11105110 / Old Material Number: 69305100D



- Suitable for standard mobile injector units and injector baskets
- Suitable e.g. for measuring cylinders, flasks and test tubes
- Including nozzle cap, Ø 5 mm

Product category	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PLW 6011	•
PLW 6111	•
Product category	
Injector nozzle	•
Product properties	
Material	Stainless steel, Plastic
Colour	Blue, Stainless Steel
Area of application	
Reprocessing of laboratory glassware	•
Dimensions and weight	
External dimensions, net height in mm	54
External dimensions, net width in mm	54
External dimensions, net depth in mm	115
External dimensions, gross height in mm	12
External dimensions, gross width in mm	80
External dimensions, gross depth in mm	120
Net weight in kg	0.01
Gross weight in kg	0.02