

APLW 052

Injector nozzle

with plastic support, star-shaped with clamp, Ø 4, length 175 mm.

EAN: 4002516117278 / Material number: 11105120 / Old Material Number: 69305200D



- Suitable for standard mobile injector units and injector baskets
- Suitable e.g. for measuring cylinders, flasks and test tubes
- Including nozzle cap, Ø 5 mm

Product classification	
Washer-disinfectors, laboratory	•
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PLW 6011	•
PLW 6111	•
PLW 8505	•
Product category	
Injector nozzle	•
Product properties	
Material	Stainless steel, Plastic
Colour	Blue, Stainless Steel
Area of application	
Reprocessing of laboratory glassware	•
Dimensions and weight	
External dimensions, net height in mm	54
External dimensions, net width in mm	54
External dimensions, net depth in mm	178
External dimensions, gross height in mm	12
External dimensions, gross width in mm	100
External dimensions, gross depth in mm	150
Net weight in kg	0.02
Gross weight in kg	0.02