

APLW 066
 Injector nozzle
 with holder, Ø 6, length 175 mm.

EAN: 4002516116325 / Material number: 11105260 / Old Material Number: 69306600D



- Suitable for standard mobile injector units and injector baskets
- Suitable e.g. for measuring cylinders, flasks and test tubes
- Including nozzle cap, Ø 10 mm

| Product classification | |
|------------------------------------------------|-----------------------------------------|
| Large-capacity washer-disinfectors, laboratory | • |
| Machine affiliation | |
| PLW 6011 | • |
| PLW 6111 | • |
| Product category | |
| Injector nozzle | • |
| Product properties | |
| Material | Wire, rilsan-coated, Stainless steel |
| Colour | Blue, Stainless Steel, White |
| Area of application | |
| Reprocessing of laboratory glassware | • |
| Dimensions and weight | |
| External dimensions, net height in mm | 75 |
| External dimensions, net width in mm | 75 |
| External dimensions, net depth in mm | 182 |
| External dimensions, gross height in mm | 55 |
| External dimensions, gross width in mm | 180 |
| External dimensions, gross depth in mm | 250 |
| Net weight in kg | 0.07 |
| Gross weight in kg | 0.08 |