

APLW 069

Injector nozzle

with stainless steel support, Ø 8, length 255 mm.

EAN: 4002516116356 / Material number: 11105290 / Old Material Number: 69306900D



- Suitable for standard mobile injector units and injector baskets
- Suitable e.g. for measuring cylinders, flasks and test tubes
- Contact surface at tip of nozzle, Ø 17 mm

Product category	
Large-chamber laboratory washers	•
Machine affiliation	
PLW 6011	•
PLW 6111	•
Product category	
Injector nozzle	•
Area of application	
Reprocessing of laboratory glassware	•
Product properties	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Grey
Dimensions and weight	
External dimensions, net height in mm	87
External dimensions, net width in mm	87
External dimensions, net depth in mm	263
External dimensions, gross height in mm	12
External dimensions, gross width in mm	250
External dimensions, gross depth in mm	350
Net weight in kg	0.08
Gross weight in kg	0.10