

APLW 070

Injector nozzle

with stainless steel support, Ø 8, length 320 mm.



- Suitable for standard mobile injector units and injector baskets
- Suitable e.g. for measuring cylinders, flasks and test tubes
- Contact surface at tip of nozzle, Ø 17 mm

EAN: 4002516116363 / Material number: 11105300 / Old Material Number: 69307000D

Product classification	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PLW 6011	•
PLW 6111	•
Product category	
Injector nozzle	•
Product properties	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Grey
Area of application	
Reprocessing of laboratory glassware	•
Dimensions and weight	
External dimensions, net height in mm	87
External dimensions, net width in mm	87
External dimensions, net depth in mm	330
External dimensions, gross height in mm	108
External dimensions, gross width in mm	117
External dimensions, gross depth in mm	305
Net weight in kg	0.10
Gross weight in kg	0.19