

APLW 071

Injector nozzle

w. holder for meas. cylinders up to Ø 105 mm & injector nozzles Ø 8, l. 320 mm.

EAN: 4002516116370 / Material number: 11105310 / Old Material Number: 69307100D



- Suitable for standard mobile injector units and injector baskets
- Suitable e.g. for measuring cylinders, flasks and test tubes
- Contact surface at tip of nozzle, Ø 17 mm

<b>Product classification</b>	
Large-capacity washer-disinfectors, laboratory	•
<b>Machine affiliation</b>	
PLW 6011	•
PLW 6111	•
<b>Product category</b>	
Injector nozzle	•
<b>Product properties</b>	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Grey
<b>Area of application</b>	
Reprocessing of laboratory glassware	•
<b>Dimensions and weight</b>	
External dimensions, net height in mm	120
External dimensions, net width in mm	120
External dimensions, net depth in mm	330
External dimensions, gross height in mm	125
External dimensions, gross width in mm	165
External dimensions, gross depth in mm	350
Net weight in kg	0.24
Gross weight in kg	0.38