

APLW 072

Blind stopper

For sealing injector positions on load carriers.



EAN: 4002516116387 / Material number: 11105320 / Old Material Number: 69307200D

Product category	
Washer-disinfectors, laboratory	•
Large-capacity washer-disinfectors, laboratory	•
Washer-disinfectors, medical	•
Machine affiliation	
PWD 8534	•
PLW 6011	•
PLW 6111	•
Product category	
Other accessories	•
Product properties	
Material	Stainless steel
Colour	Stainless Steel
Area of application	
Reprocessing of laboratory glassware	•
Dimensions and weight	
External dimensions, net height in mm	10
External dimensions, net width in mm	10
External dimensions, net depth in mm	15
External dimensions, gross height in mm	12
External dimensions, gross width in mm	80
External dimensions, gross depth in mm	120
Net weight in kg	0.00
Gross weight in kg	0.00