


APLW 083
Support
for injector nozzle, height 130 mm.



EAN: 4002516116592 / Material number: 11105430 / Old Material Number: 69308300D

Product classification	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PLW 6011	•
PLW 6111	•
Product category	
Other accessories	•
Product properties	
Material	Stainless steel, Plastic
Colour	Blue, Stainless Steel
Area of application	
Reprocessing of laboratory glassware	•
Dimensions and weight	
External dimensions, net height in mm	50
External dimensions, net width in mm	50
External dimensions, net depth in mm	133
External dimensions, gross height in mm	55
External dimensions, gross width in mm	180
External dimensions, gross depth in mm	250
Net weight in kg	0.02
Gross weight in kg	0.03