

APLW 085

Bottle neck holder

Ø 33 mm, for injector nozzles Ø 6 mm.



**Image**  
in progress

EAN: 4002516116615 / Material number: 11105450 / Old Material Number: 69308500D

<b>Product category</b>	
Large-capacity washer-disinfectors, laboratory	•
<b>Machine affiliation</b>	
PLW 6011	•
PLW 6111	•
<b>Product category</b>	
Other accessories	•
<b>Product properties</b>	
Material	Plastic
Colour	White
<b>Area of application</b>	
Reprocessing of laboratory glassware	•
<b>Dimensions and weight</b>	
External dimensions, net height in mm	25
External dimensions, net width in mm	25
External dimensions, net depth in mm	42
External dimensions, gross height in mm	12
External dimensions, gross width in mm	120
External dimensions, gross depth in mm	170
Net weight in kg	0.02
Gross weight in kg	0.02