

APLW 093
 Irrigation sleeves
 with cap for pipettes, maximum Ø 11 mm



EAN: 4002516116721 / Material number: 11105530 / Old Material Number: 69309300D

Product category	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PLW 6011	•
PLW 6111	•
Product category	
Holder	•
Other accessories	•
Product properties	
Material	Stainless steel, Silicone
Colour	Stainless Steel, White
Area of application	
Reprocessing of laboratory glassware	•
Dimensions and weight	
External dimensions, net height in mm	25
External dimensions, net width in mm	25
External dimensions, net depth in mm	118
External dimensions, gross height in mm	12
External dimensions, gross width in mm	80
External dimensions, gross depth in mm	120
Net weight in kg	0.05
Gross weight in kg	0.05