


A 301/5

Module

for optimum loading of laboratory glassware – can be equipped with 18 nozzles.



- Suitable for lab. glassware incl. Erlenmeyer and round flasks.
- Equipped with 9 A 842 nozzles and 9 A 843 nozzles
- Part of the EasyLoad system
- Nozzle spacing 78 mm across the width, 78 mm across the depth

EAN: 4002516188971 / Material number: 11116950 / Old Material Number: 69630105D

<b>Product category</b>	
Washer-disinfectors, laboratory	•
Large-capacity washer-disinfectors, laboratory	•
<b>Machine affiliation</b>	
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
<b>Product category</b>	
Insert/module for narrow-necked glassware	•
Insert for various flasks	•
<b>Product properties</b>	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Black
Part of the EasyLoad system	•
<b>Area of application</b>	
Reprocessing of laboratory glassware	•
<b>Capacity</b>	
Narrow-necked glassware [number]	18
<b>Dimensions and weight</b>	
External dimensions, net height in mm	203
External dimensions, net width in mm	232
External dimensions, net depth in mm	480
External dimensions, gross height in mm	200
External dimensions, gross width in mm	245
External dimensions, gross depth in mm	515
Net weight in kg	2.16
Gross weight in kg	3.17
<b>Companion pack</b>	
A 842, injector nozzle Ø 4 mm, length 90 mm [number]	9
A 843, injector nozzle Ø 4 mm, length 185 mm [number]	9