

APWD 072

Injector nozzle

for optimum reprocessing of lumened instruments, particularly suction cups.

EAN: 4002516122289 / Material number: 11127460 / Old Material Number: 68407200D



- Suitable for the APWD 061 basket
- Using adapter APWD 073
- Figure shows APWD 072 with APWD 073 adapter

Product classification	
Thermal disinfectors	•
Machine affiliation	
PWD 8531	•
PWD 8532	•
Product category	
Injector nozzle	•
Product properties	
Material	Stainless steel
Colour	Stainless Steel
Area of application	
Reprocessing of dental instruments	•
Reprocessing of lumened instruments	•
Dimensions and weight	
External dimensions, gross height in mm	12
External dimensions, gross width in mm	80
External dimensions, gross depth in mm	120
Net weight in kg	0.00
Gross weight in kg	0.00