

APWD 073

Adapter

for additional positioning of APWD 072 injector nozzles in the APWD 061 basket.

EAN: 4002516122296 / Material number: 11127500 / Old Material Number: 68407300D



Product classification	
Thermal disinfectors	•
Machine affiliation	
PWD 8531	•
PWD 8532	•
Product category	
Adapter	•
Other accessories	•
Product properties	
Material	Stainless steel
Colour	Stainless Steel
Area of application	
Reprocessing of dental instruments	•
Reprocessing of lumened instruments	•
Dimensions and weight	
External dimensions, gross height in mm	12
External dimensions, gross width in mm	80
External dimensions, gross depth in mm	120
Net weight in kg	0.00
Gross weight in kg	0.01