


A 851
Holder rack
For injector module A 621.



- Arrangement: 3 x 7 fields
- For use on A 621 injector module
- For the optimum positioning of laboratory glassware
- Part of the EasyLoad system

EAN: 4002516201915 / Material number: 11236640 / Old Material Number: 69685100D

Product category	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PLW 8615	•
PLW 8616	•
PLW 8617	•
Product category	
Other inserts/other modules	•
Product properties	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Black
Number of compartments	21
Compartment width in mm	74
Compartment length in mm	77
Part of the EasyLoad system	•
Area of application	
Reprocessing of laboratory glassware	•
Dimensions and weight	
External dimensions, net height in mm	119
External dimensions, net width in mm	243
External dimensions, net depth in mm	580
External dimensions, gross height in mm	135
External dimensions, gross width in mm	285
External dimensions, gross depth in mm	586
Net weight in kg	0.67
Gross weight in kg	1.07