


A 852  
Holder rack  
For injector module A 622.



- Arrangement: 4 x 9 fields
- For use on A 622 injector module
- For the optimum positioning of laboratory glassware
- Part of the EasyLoad system

EAN: 4002516201922 / Material number: 11236650 / Old Material Number: 69685200D

<b>Product category</b>	
Large-capacity washer-disinfectors, laboratory	•
<b>Machine affiliation</b>	
PLW 8615	•
PLW 8616	•
PLW 8617	•
<b>Product category</b>	
Other inserts/other modules	•
<b>Product properties</b>	
Material	Wire, rilsan-coated, Stainless steel, Plastic
Colour	Stainless Steel, Grey, Black
Number of compartments	36
Compartment width in mm	59
Compartment length in mm	58
Part of the EasyLoad system	•
<b>Area of application</b>	
Reprocessing of laboratory glassware	•
<b>Dimensions and weight</b>	
External dimensions, net height in mm	76
External dimensions, net width in mm	256
External dimensions, net depth in mm	562
External dimensions, gross height in mm	75
External dimensions, gross width in mm	330
External dimensions, gross depth in mm	680
Net weight in kg	0.61
Gross weight in kg	0.99