

A 12/2

Insert

For the gentle cleaning of delicate load items by interrupting the water jet.



- With perforations 7 x 7 mm
- For use in upper and lower basket
- Covers half of basket
- Not for use in PG 8536

EAN: 4002516237068 / Material number: 11238030 / Old Material Number: 69501205D

Product category	
Washer-disinfectors, laboratory	•
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PG 8535	•
PG 8536	•
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
PLW 8505	•
Product category	
Underlay	•
Retaining frame	•
Product properties	
Material	Stainless steel
Colour	Stainless Steel
Area of application	
Reprocessing of laboratory glassware	•
Dimensions and weight	
External dimensions, net height in mm	237
External dimensions, net width in mm	448
External dimensions, net depth in mm	14
External dimensions, gross height in mm	65
External dimensions, gross width in mm	255
External dimensions, gross depth in mm	475
Net weight in kg	0.39
Gross weight in kg	0.57