


A 11/2

Insert

For the gentle cleaning of delicate load items by interrupting the water jet.



- With perforations 7 x 7 mm
- For use in upper and lower basket

EAN: 4002516340119 / Material number: 11239130 / Old Material Number: 69501103D

| Product classification | |
|--|-----------------|
| Washer-disinfectors, laboratory | • |
| Large-capacity washer-disinfectors, laboratory | • |
| Machine affiliation | |
| PG 8504 | • |
| PG 8583 | • |
| PG 8583 CD | • |
| PG 8593 | • |
| PLW 8505 | • |
| Product category | |
| Underlay | • |
| Retaining frame | • |
| Product properties | |
| Material | Stainless steel |
| Colour | Stainless Steel |
| Area of application | |
| Reprocessing of laboratory glassware | • |
| Dimensions and weight | |
| External dimensions, net height in mm | 1 |
| External dimensions, net width in mm | 476 |
| External dimensions, net depth in mm | 443 |
| External dimensions, gross height in mm | 70 |
| External dimensions, gross width in mm | 505 |
| External dimensions, gross depth in mm | 430 |
| Net weight in kg | 0.71 |
| Gross weight in kg | 1.08 |