

A 320

Insert

For the gentle cleaning of delicate load items by interrupting the water jet.



- With perforations 7 x 7 mm

EAN: 4002516237082 / Material number: 11295510 / Old Material Number: 69632000D

| | |
|--|-----------------|
| Product category | |
| Large-capacity washer-disinfectors, laboratory | • |
| Machine affiliation | |
| PLW 8615 | • |
| PLW 8616 | • |
| PLW 8617 | • |
| Product category | |
| Retaining frame | • |
| Product properties | |
| Material | Stainless steel |
| Colour | Stainless Steel |
| Area of application | |
| Reprocessing of laboratory glassware | • |
| Dimensions and weight | |
| External dimensions, net height in mm | 10 |
| External dimensions, net width in mm | 237 |
| External dimensions, net depth in mm | 536 |
| External dimensions, gross height in mm | 60 |
| External dimensions, gross width in mm | 245 |
| External dimensions, gross depth in mm | 560 |
| Net weight in kg | 0.47 |
| Gross weight in kg | 0.63 |