

A 320

Insert

For the gentle cleaning of delicate load items by interrupting the water jet.



- With perforations 7 x 7 mm

EAN: 4002516237082 / Material number: 11295510 / Old Material Number: 69632000D

Product category	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PLW 8615	•
PLW 8616	•
PLW 8617	•
PLW 7111	•
Product category	
Retaining frame	•
Product properties	
Material	Stainless steel
Colour	Stainless Steel
Area of application	
Reprocessing of laboratory glassware	•
Dimensions and weight	
External dimensions, net height in mm	10
External dimensions, net width in mm	237
External dimensions, net depth in mm	536
External dimensions, gross height in mm	60
External dimensions, gross width in mm	245
External dimensions, gross depth in mm	560
Net weight in kg	0.47
Gross weight in kg	0.63