

DOS K 85/1 flex

Dispenser module for 5 litre canister

with level indicator for the automatic dispensing of liquid media.



- Dispenser monitoring – high process security as per EN ISO 15883
- For dispensing liquid cleaning agent/neutralising agent
- 214 mm siphon for 5 l media canisters with flexible cap
- Siphon conversion kit for 10–30 l media canisters available

EAN: 4002516365297 / Material number: 11498920 / Old Material Number: 69747053D

| | |
|---|---------|
| Product classification | |
| Washer-disinfectors, laboratory | • |
| Thermal disinfectors | • |
| Washer-disinfectors, medical | • |
| Machine affiliation | |
| PG 8581 | • |
| PG 8562 | • |
| PG 8582 | • |
| PG 8583 | • |
| PG 8591 | • |
| PG 8592 | • |
| PG 8593 | • |
| Product category | |
| Dispenser system | • |
| Product properties | |
| Production rate in ml/60 s | 55 |
| Length of suction lance in mm | 214 |
| Standard electrical connection | |
| Number of phases | 1 |
| Neutral conductor | N |
| AC/DC | AC |
| Voltage in V | 200-230 |
| Frequency in Hz | 50 |
| Total rated load in kW | 0.01 |
| Fuse rating in A | 0 |
| Dimensions and weight | |
| External dimensions, net height in mm | 0 |
| External dimensions, net width in mm | 0 |
| External dimensions, net depth in mm | 0 |
| External dimensions, gross height in mm | 160 |
| External dimensions, gross width in mm | 198 |
| External dimensions, gross depth in mm | 478 |
| Net weight in kg | 1.52 |
| Gross weight in kg | 1.74 |
| Standards, certs, symbols | |
| CE | • |