

A 830/1

Insert

ideal holder for separable mouth specula.



- For use in insert E 337/1 and E 802/1
- 3 compart. Ø 8 mm, 14 compart. 10 x 19 mm, 2 compart. 8 x 19 mm
- Mouth specula handles can be placed in compartments of inserts
- Suitable for mouth specula  $\geq 3$

EAN: 4002516429203 / Material number: 11729480 / Old Material Number: 68683001D

<b>Product category</b>	
Thermal disinfectors	•
<b>Machine affiliation</b>	
PG 8528	•
G 7831	•
PG 8536	•
PG 8581	•
PG 8591	•
PWD 8531	•
PWD 8532	•
<b>Product category</b>	
Insert/module for dental instruments	•
<b>Product properties</b>	
Material	Plastic
Colour	Black
Number of compartments	16
<b>Area of application</b>	
Reprocessing of dental instruments	•
<b>Dimensions and weight</b>	
External dimensions, net height in mm	14
External dimensions, net width in mm	49
External dimensions, net depth in mm	159
External dimensions, gross height in mm	67
External dimensions, gross width in mm	170
External dimensions, gross depth in mm	230
Net weight in kg	0.02
Gross weight in kg	0.02