

A 813/1

Adapter

optimum reprocessing of transmission instruments with external spray channel.



- For use in A 105/1 upper baskets and in APWD 061 baskets
- Two-piece, comprising an upper and lower part incl. accessories

EAN: 4002516432555 / Material number: 11763880 / Old Material Number: 68681301D

Product classification	
Thermal disinfectors	•
Machine affiliation	
PG 8581	•
PG 8591	•
PWD 8531	•
PWD 8532	•
Product category	
Holder and/or adapter for transmission instruments	•
Product properties	
Material	Stainless steel
Colour	Stainless Steel
Area of application	
Reprocessing of dental instruments	•
Capacity	
Transmission instruments [number]	4
Dimensions and weight	
External dimensions, net height in mm	120
External dimensions, net width in mm	120
External dimensions, net depth in mm	300
External dimensions, gross height in mm	108
External dimensions, gross width in mm	117
External dimensions, gross depth in mm	305
Net weight in kg	0.19
Gross weight in kg	1.00
Companion pack	
E 362 blind stopper [number]	3
Nozzle with silicone hose [number]	4