

WA UC 1201 LB

Miele UltraColor refill bag, 1.265 l
for coloureds and dark garments.

EAN: 4002516430919 / Material number: 11764630 / Old Material Number: 11997207EU1



- Excellent spot removal even at low temperatures
- Excellent spot removal even at low temperatures
- With colour-protection formula for brilliant colours that last
- Free from boron compounds, phosphates and micro-plastics
- Refill bag with practical dispenser
- Highly efficient – for 23 wash loads
- Simply perfect laundry care. With Miele.