

## XKM 3200 WL LMD

### Communication module

for connection to Miele MOVE, the Segosoft Miele Edition or other software.



- State-of-the-art connection to the Miele Cloud Service
- Secure transmission over encrypted WiFi or Ethernet
- REST interface with data provision in JSON format
- Flexible operation on local networks
- Up to 30 process records are stored in the appliance
- This module is not compatible with DataDiary

EAN: 4002516567226 / Material number: 11795280 / Old Material Number: 68822103D

Product classification	
Washer-disinfectors, laboratory	•
Thermal disinfectors	•
Washer-disinfectors, medical	•
Machine affiliation	
PG 8581	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8583	•
PG 8583 CD	•
PG 8591	•
PG 8592	•
PG 8593	•
Product category	
Process documentation/operating documentation/communication	•
Product properties	
Material	Plastic
Colour	Grey
Area of application	
Reprocessing of pipettes/phials	•
Reprocessing of laboratory glassware	•
Reprocessing of anaesthetic instruments	•
Reprocessing of MIS instruments	•
Reprocessing of ophthalmology instruments	•
Reprocessing of robotic instruments	•
Reprocessing of GYN instruments	•
Reprocessing of ENT instruments	•
Reprocessing of dental instruments	•
Reprocessing of theatre shoes	•
Reprocessing of basins	•
Reprocessing of kidney dishes	•
Reprocessing of baby bottles	•
Reprocessing of surgical instruments	•
Reprocessing of lumened instruments	•
Reprocessing of orthopaedic surgery instruments	•
Reprocessing of intramedullary nails/bone marrow drills	•
Reprocessing of trays	•
Reprocessing of mesh trays	•
Reprocessing of different utensils	•
Dimensions and weight	
External dimensions, net height in mm	39
External dimensions, net width in mm	88
External dimensions, net depth in mm	110
External dimensions, gross height in mm	140
External dimensions, gross width in mm	190
External dimensions, gross depth in mm	295
Net weight in kg	0.61
Gross weight in kg	0.86

XKM 3200 WL LMD

Communication module

for connection to Miele MOVE, the Segosoft Miele Edition or other software.

EAN: 4002516567226 / Material number: 11795280 / Old Material Number: 68822103D



Standards, certs, symbols	
CE	•
FCC	•
UKCA	•
EAC	•