


APLW 152  
Basket  
for various inserts.



- For use on level 1
- Vertical clearance dependent on other baskets used
- Clearances: width 507 mm, depth 545 mm

EAN: 4002516526001 / Material number: 11854310 / Old Material Number: 69615200D

<b>Product category</b>	
Large-capacity washer-disinfectors, laboratory	•
<b>Machine affiliation</b>	
PLW 7111	•
PLW 7111 [Safety]	•
<b>Product category</b>	
Carrier/basket carrier	•
<b>Product properties</b>	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Black
<b>Dimensions and weight</b>	
External dimensions, net height in mm	180
External dimensions, net width in mm	559
External dimensions, net depth in mm	568
External dimensions, gross height in mm	180
External dimensions, gross width in mm	570
External dimensions, gross depth in mm	625
Net weight in kg	2.33
Gross weight in kg	3.41