

WA UC 1501 L

UltraColor liquid detergent 1.5 l

highly effective liquid detergent for colours and dark / black items.

EAN: 4002516531500 / Material number: 11979350 / Old Material Number: 11997231EU2



- Excellent spot removal even at low temperatures
- Contains 7 enzymes for ultimate cleanliness
- Bottles made from 100 % recycled plastic – eco-friendly
- ColorProtect for long-lasting colour intensity
- Free from microplastics and biodegradable
- Long-lasting freshness with the Miele Aqua fragrance