

FA DF 452 L  
Set of 3x Miele DryFresh™  
for 150 drying cycles – Freshen-up and fragrance.

EAN: 4002516561200 / Material number: 12026040 / Old Material Number: 11997236EU1



- With Freshplex™ to neutralise unpleasant odours
- Freshens up dry laundry in the “Timed drying” programme
- Beautifully fragrant laundry – without fabric conditioner
- Can be combined with other Miele fragrance flacons
- Specially developed for Miele by perfumers from in the Provence
- For all Miele T1 tumble dryers with FragranceDos