

APWD 328

Insert

expansion set for reprocessing dental nitrous oxide systems.

EAN: 4002516611028 / Material number: 12036450 / Old Material Number: 68632800D



- Capacity expansion to up to a total of 5 double hose sets
- 5x clip and swivel plates

Product classification	
Thermal disinfectors	•
Machine affiliation	
PG 8581	•
PG 8591	•
Product category	
Other inserts	•
Product properties	
Material	Stainless steel
Colour	Stainless Steel
Area of application	
Reprocessing of dental instruments	•
Reprocessing of different utensils	•
Dimensions and weight	
External dimensions, net height in mm	77
External dimensions, net width in mm	90
External dimensions, net depth in mm	100
External dimensions, gross height in mm	80
External dimensions, gross width in mm	95
External dimensions, gross depth in mm	110
Net weight in kg	0.26
Gross weight in kg	0.32