

Complete C3 Cat & Dog  
Cylinder vacuum cleaner  
With turbobrush and odour filter; ideal for pet lovers.



- Reliably removes hair and lint with its turbo brush.
- Odours are neutralised more effectively – Active AirClean filter
- Easy cleaning of larger areas thanks to 12 m operating radius
- Power Efficiency Motor with improved efficiency
- Convenient and versatile with three integrated accessories

EAN: 4002516599401 / Material number: 12137590 / Old Material Number: 41GEF532GB

<b>Construction type</b>	
Cylinder vacuum cleaner with dustbag	•
Product brand	Complete C3
Motor type	PowerLine
Product type	SGEF5
Appliance Colour	Mango red
<b>User convenience</b>	
Electronic suction power control	+/- footswitches
Illuminated power indicator	•
Ergonomic carrying handle	•
Ergonomic handle	•
Silence system	•
Suction tube	Comfort telescopic tube
Park system for pauses during vacuuming	•
Park system	on both sides
Comfort cable rewind	•
<b>Efficiency and sustainability</b>	
Nominal power consumption in W	890
<b>Mobility</b>	
Castors with steel shaft, gentle on floors	•
3-point manoeuvrability	•
Operating radius in m	12
Cable length in m	8.5
<b>Filter system</b>	
Vacuum cleaner bag type	HyClean GN
Dustbag volume in l	4.5
Dustbag change indicator	•
Exhaust filter	Active AirClean filter
Motor protection filter	•
<b>Safety</b>	
Locking system for floorhead and handle	•
Automatic dustbag positioning	•
All-round furniture protection strip	•
<b>Technical data</b>	
Max. power	890
Volume in dB	76.0
Width in mm	496
Height in mm	227
Weight in kg	7.3
<b>Accessories included</b>	
Three-piece accessory set	Integrated
Upholstery nozzle	•
Crevice nozzle	•
Dusting brush with synthetic bristles	•
SBD 365-3 floorhead	•
Turbobrush	STB 305-3
Turbobrush STB 305-3	•
Filter	Active AirClean filter
Cleaning tool	No