

APLW 631

Insert

for the optimum loading of breathing masks.



- Safe reprocessing of up to 4 breathing masks
- With 4 small-part baskets

EAN: 4002516726036 / Material number: 12230130 / Old Material Number: 69663100D

Product category	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PLW 8615 [Safety]	•
PLW 8616 [Safety]	•
PLW 8617 [Safety]	•
PLW 7111 [Safety]	•
Product category	
Insert for breathing masks	•
Product properties	
Material	Stainless steel
Colour	Stainless Steel
Area of application	
Reprocessing of breathing masks	•
Capacity	
Breathing masks [number]	4
Dimensions and weight	
External dimensions, gross height in mm	190
External dimensions, gross width in mm	265
External dimensions, gross depth in mm	550
Net weight in kg	2.99
Gross weight in kg	2.99