

APLW 873

Hangers

with carabiner hooks for colour coding.



- Set of 6 hangers in different colours and 6 carabiner hooks
- For use with APFD 200, APFD 201, APLW 630 and APLW 631

EAN: 4002516726050 / Material number: 12230230 / Old Material Number: 69687300D

|  |  |
|--|--|
| <b>Product category</b>                        |  |
| Washer-disinfectors, laboratory                | •  |
| Large-capacity washer-disinfectors, laboratory | •  |
| <b>Machine affiliation</b>                     |  |
| PG 8063  | •  |
| PG 8593  | •  |
| PLW 8615 [Safety]                              | •  |
| PLW 8616 [Safety]                              | •  |
| PLW 8617 [Safety]                              | •  |
| PLW 7111 [Safety]                              | •  |
| <b>Product category</b>                        |  |
| Other accessories                              | •  |
| <b>Product properties</b>                      |  |
| Material                                       | Stainless steel, Plastic                                   |
| Colour   | Blue, Stainless Steel, Yellow, Green, Black, White, Purple |
| <b>Area of application</b>                     |  |
| Reprocessing of breathing air regulators       | •  |
| Reprocessing of breathing masks                | •  |
| <b>Dimensions and weight</b>                   |  |
| External dimensions, gross height in mm        | 12   |
| External dimensions, gross width in mm         | 80   |
| External dimensions, gross depth in mm         | 120  |
| Net weight in kg                               | 0.03   |
| Gross weight in kg                             | 0.03   |