


APFD 421

Cutlery basket

for additional cutlery capacity in the lower basket of ProfiLine dishwashers.

EAN: 4002516678687 / Material number: 12256760 / Old Material Number: 67042101D



- Easily retrofittable, in addition to the cutlery tray
- Particularly suitable for bulky items such as whisks
- Cutlery basket can be positioned conveniently in the lower basket

<b>Product category</b>	
Fresh-water dishwashers	•
<b>Machine affiliation</b>	
PFD 100 (SmartBiz)	•
PFD 101	•
PFD 101 i	•
PFD 101 U	•
PFD 101 U/i	•
PFD 102 i	•
PFD 103 SCi XXL	•
PFD 104 SCVi XXL	•
<b>Product category</b>	
Insert for cutlery	•
<b>Product properties</b>	
Material	Plastic
Colour	White
<b>Dimensions and weight</b>	
External dimensions, gross height in mm	240
External dimensions, gross width in mm	190
External dimensions, gross depth in mm	280
Net weight in kg	0.47
Gross weight in kg	0.80