


APFD 222

Basket

for optimum loading of reusable tumblers.



- With locking mechanism to secure cups during washing
- For 25 medication cups with a maximum height of 200 mm
- For cups with a diameter of min. 41–60 mm and max. 90 mm

EAN: 4002516681595 / Material number: 12294760 / Old Material Number: 67022201D

<b>Product category</b>	
Tank dishwashers	•
Throughfeed tank dishwashers	•
<b>Machine affiliation</b>	
PTD 901	•
PFD 400	•
PFD 400 U	•
PFD 401	•
PFD 401 U	•
PFD 402	•
PFD 402 U	•
PFD 404	•
PFD 404 U	•
PFD 405	•
PFD 405 U	•
PFD 407	•
PFD 407 U	•
PFD 408 U	•
PTD 702	•
PTD 703	•
PTD 704	•
<b>Product category</b>	
Basket for beakers	•
<b>Capacity</b>	
Beakers [number]	25
<b>Dimensions and weight</b>	
External dimensions, net height in mm	214
External dimensions, net width in mm	504
External dimensions, net depth in mm	514
External dimensions, gross height in mm	223
External dimensions, gross width in mm	563
External dimensions, gross depth in mm	563
Net weight in kg	2.86
Gross weight in kg	3.85