

APLW 880

Retrofitting kit

for Rec in, fourth water inlet, including booster pump for PLW 8636.



- Required in order to pump a medium out of a buffer container
- For use with an existing on-site recycling tank
- Can be used with existing on-site container, e.g. IBC tank
- Installation is carried out by Customer Service, additional costs possible

EAN: 4002516745532 / Material number: 12426490 /
Old Material Number: 69688000D

Product category

- | | |
|---|---|
| Washer-disinfectors, laboratory | • |
| Cleaning machines for industrial parts cleaning | • |

Machine affiliation

- | | |
|--------------------------|---|
| PLW 8636 [IPC] | • |
| PLW 8636 [LAB / LAB MON] | • |

Product category

- | | |
|----------------|---|
| Conversion kit | • |
|----------------|---|

Dimensions and weight

- | | |
|---|------|
| External dimensions, gross height in mm | 123 |
| External dimensions, gross width in mm | 315 |
| External dimensions, gross depth in mm | 505 |
| Net weight in kg | 4.36 |
| Gross weight in kg | 4.84 |