

APLW 881

Retrofitting kit

for Rec out, additional water drain, incl. booster pump for PLW 8636.



- Required in order to pump a medium into a buffer container
- For use with an existing on-site recycling tank
- Can be used with existing on-site container, e.g. IBC tank
- Installed by Customer Service

EAN: 4002516738978 / Material number: 12426550 / Old Material Number: 69688100D

Product category

Washer-disinfectors, laboratory

•

Cleaning machines for industrial parts cleaning

•

Machine affiliation

PLW 8636 [IPC]

•

PLW 8636 [LAB / LAB MON]

•

Product category

Conversion kit

•

Dimensions and weight

External dimensions, gross height in mm

123

External dimensions, gross width in mm

315

External dimensions, gross depth in mm

505

Net weight in kg

4.14

Gross weight in kg

4.62