

APLW 885

Retrofitting kit

extra disp. membrane pump & suction lance, ext. 10–20 l canisters for PLW 8636.

EAN: 4002516745525 / Material number: 12426770 / Old Material Number: 69688500D



Product category	
Washer-disinfectors, laboratory	•
Cleaning machines for industrial parts cleaning	•
Machine affiliation	
PLW 8636 [IPC]	•
PLW 8636 [LAB / LAB MON]	•
Product category	
Conversion kit	•
Dispenser system	•

- For 10–20 l canisters, positioned outside the machine
- The suction lance features level detection
- Installed by Customer Service