


APLW 886

Retrofitting kit

for dispensing monitoring by oval-wheel counter, for installation in PLW 8636.

EAN: 4002516745433 / Material number: 12426910 / Old Material Number: 69688600D



- More precise quantity measurement using oval-wheel counters
- Installation is carried out by Customer Service, additional costs possible

Product category

Washer-disinfectors, laboratory

•

Cleaning machines for industrial parts cleaning

•

Machine affiliation

PLW 8636 [IPC]

•

PLW 8636 [LAB / LAB MON]

•

Product category

Conversion kit

•

Dispenser system

•

Dimensions and weight

External dimensions, gross height in mm

80

External dimensions, gross width in mm

95

External dimensions, gross depth in mm

110

Net weight in kg

0.03

Gross weight in kg

0.08