

APLW 887

Retrofitting kit

additional suction lance, for external 10–20 l canisters for PLW 8636.

EAN: 4002516745082 / Material number: 12426920 / Old Material Number: 69688700D



Product category	
Washer-disinfectors, laboratory	•
Cleaning machines for industrial parts cleaning	•
Machine affiliation	
PLW 8636 [IPC]	•
PLW 8636 [LAB / LAB MON]	•
Product category	
Conversion kit	•
Dispenser system	•

- For 10–20 l canisters, positioned outside the machine
- The suction lance features level detection
- Installation is carried out by Customer Service, additional costs possible