

## A 315/1 Module

For increased capacity and reprocessing of lumened instruments.



- Increases the injector capacity from 22 to 44
- Can be used in all A 151 lower baskets available on the market
- Compatible with all standard dental & ENT accessories
- Injector nozzles not included
- A 800 tubular filter for use in either injector manifold (not included)
- Clearances when used in A 151 lower basket whilst at the same time using A 105/1

EAN: 4002516881186 / Material number: 12531630 / Old Material Number: 69631501D

Product category	
Thermal disinfectors	•
Washer-disinfectors, medical	•
Machine affiliation	
PG 8581	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8591	•
PG 8592	•
PWD 8682	•
PWD 8682 C	•
PWD 8682 CD	•
PWD 8692	•
Product category	
Insert for ENT instruments	•
Insert/module for dental instruments	•
Other inserts/other modules	•
Product properties	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Black
Option to use a central filter	•
Area of application	
Reprocessing of ENT instruments	•
Reprocessing of dental instruments	•
Reprocessing of lumened instruments	•
Capacity	
ENT instruments [number]	22
Transmission instruments [number]	22
Siphons [number]	22
Lumened instruments [number]	22
Dimensions and weight	
External dimensions, gross height in mm	165
External dimensions, gross width in mm	505
External dimensions, gross depth in mm	550
Net weight in kg	2.90
Gross weight in kg	3.66
Companion pack	
A 838 set of spanners [number]	1