


APWD 153

Basket

For holding standard DIN mesh trays as well as various inserts.



EAN: 4002516847540 / Material number: 12537450 / Old Material Number: 69615300D

Product category	
Thermal disinfectors	•
Washer-disinfectors, medical	•
Machine affiliation	
PWD 8682 C	•
Product category	
Carrier/basket carrier	•
Product properties	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Black
Area of application	
Reprocessing of GYN instruments	•
Reprocessing of ENT instruments	•
Reprocessing of dental instruments	•
Reprocessing of VET instruments	•
Reprocessing of podiatry instruments	•
Reprocessing of basins	•
Reprocessing of kidney dishes	•
Reprocessing of surgical instruments	•
Reprocessing of lumened instruments	•
Reprocessing of trays	•
Reprocessing of mesh trays	•
Reprocessing of small parts	•
Reprocessing of different utensils	•
Capacity	
DIN mesh trays [number]	2
Dimensions and weight	
External dimensions, net height in mm	88
External dimensions, net width in mm	530
External dimensions, net depth in mm	522
External dimensions, gross height in mm	115
External dimensions, gross width in mm	540
External dimensions, gross depth in mm	553
Net weight in kg	1.86
Gross weight in kg	2.46