

TU HyClean Pure

HyClean Pure TU vacuum cleaner bag

Miele AirClean system delivers the best in filtration performance at 99.99 %.

EAN: 4002516829751 / Material number: 12557060 / Old Material Number: 41996682INT



- 80 % recycled material and ultimate performance – unrivalled on the market
- Airflow guide and 3D expansion for long service life
- HyClose automatic self-sealing collar for even greater hygiene
- Highest level of efficiency when vacuuming with a low wattage
- Contents: 4 dustbags, 1 motor protection and 1 exhaust filter
- 16 dustbags, 4 motor protection and 4 exhaust filters