


APFD 201

Use

for the optimum loading of breathing masks.



- Suitable for holding 6 breathing masks securely
- Accessories: 1 small-part basket with 6 coloured compartments
- With separate lid for particularly delicate parts (6 compartments)

EAN: 4002516855248 / Material number: 12615430 / Old Material Number: 67020102D

| Product category | |
|---------------------------------|---|
| Fresh-water dishwashers | • |
| Machine affiliation | |
| PG 8063 | • |
| PLW 8683 [Safety] | • |
| PLW 8693 [Safety] | • |
| Product category | |
| Insert for breathing masks | • |
| Area of application | |
| Reprocessing of breathing masks | • |
| Capacity | |
| Breathing masks [number] | 6 |