


A 102

Upper basket

Height adjustable, for the optimum loading of inserts.




- Vertical clearance 205 +/- 30 mm, width 475 mm, depth 443 mm
- With integrated spray arm

EAN: 4002516855262 / Material number: 12645600 / Old Material Number: 69610201D

Product category	
Washer-disinfectors, laboratory	•
Thermal disinfectors	•
Washer-disinfectors, medical	•
Machine affiliation	
PLW 8636 [LAB / LAB MON]	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8592	•
PG 8593	•
PLW 8604	•
PLW 8683	•
PLW 8683 CD	•
PLW 8693	•
PLW 8683 [Safety]	•
PLW 8693 [Safety]	•
PWD 8682	•
PWD 8682 CD	•
PWD 8692	•
Product category	
Carrier/basket carrier	•
Upper basket	•
Product properties	
Material	Stainless steel
Colour	Stainless Steel
Area of application	
Reprocessing of laboratory glassware	•
Reprocessing of GYN instruments	•
Reprocessing of ENT instruments	•
Reprocessing of dental instruments	•
Reprocessing of theatre shoes	•
Reprocessing of basins	•
Reprocessing of kidney dishes	•
Reprocessing of surgical instruments	•
Capacity	
DIN mesh trays [number]	2
SPRI mesh trays [number]	1
Dimensions and weight	
External dimensions, net height in mm	206
External dimensions, net width in mm	528
External dimensions, net depth in mm	527
External dimensions, gross height in mm	230
External dimensions, gross width in mm	540
External dimensions, gross depth in mm	545

A 102
Upper basket
Height adjustable, for the optimum loading of inserts.



EAN: 4002516855262 / Material number: 12645600 / Old Material Number: 69610201D

Net weight in kg	3.04
Gross weight in kg	3.78