

Guard L1 Cat & Dog  
Cylinder vacuum cleaner  
powerful | attachment for removing pet hair | ideal for pet owners



- Active AirClean filter neutralises unpleasant odours
- Durable: tested for the equivalent of 20 years\*\* use
- ComfortFit: developed for the highest hygiene standards
- Individually selectable height and ergonomic Comfort handle
- For effortless cleaning: 2-in-1 upholstery and crevice nozzle
- Ideal power for every cleaning task: 4 selectable modes
- Premium design concept with stylish applications
- Turbobrush: removes stubborn dirt and animal hair

\* Compared to SBB300-3 Parquet Twister parquet brush

EAN: 4002516833215 / Material number: 12652240 /  
Old Material Number: 41UZF034GB

|  |                        |
|--|------------------------|
| <b>Construction type and design</b>        |                        |
| Appliance colour                           | Obsidian black         |
| <b>Construction type</b>                   |                        |
| Cylinder vacuum cleaner with dustbag       | •                      |
| Product brand                              | Guard L1               |
| Motor type                                 | PowerLine              |
| Product type                               | SUZFO                  |
| <b>User convenience</b>                    |                        |
| Electronic suction power control           | +/- footswitches       |
| Illuminated power indicator                | •                      |
| Suction tube                               | Suction tube X(HES)    |
| Park system for pauses during vacuuming    | •                      |
| Comfort cable rewind                       | •                      |
| Comfort handle 2.0                         | •                      |
| <b>Mobility</b>                            |                        |
| Castors with steel shaft, gentle on floors | •                      |
| 3-point manoeuvrability                    | •                      |
| Operating radius in m                      | 12                     |
| Cable length in m                          | 8.5                    |
| <b>Filter system</b>                       |                        |
| Vacuum cleaner bag type                    | HyClean Pure TU        |
| Dustbag volume in l                        | 3.5                    |
| Dustbag change indicator                   | •                      |
| Exhaust filter                             | Active AirClean filter |
| Motor protection filter                    | •                      |
| <b>Safety</b>                              |                        |
| Locking system for floorhead and handle    | •                      |
| Automatic dustbag positioning              | •                      |
| All-round furniture protection strip       | •                      |
| <b>Technical data</b>                      |                        |
| Power                                      | 890                    |
| Volume in dB                               | 74.0                   |
| Width in mm                                | 500                    |
| Height in mm                               | 239                    |
| Weight in kg                               | 6.6                    |
| <b>Accessories included</b>                |                        |
| Three-piece accessory set                  | Integrated             |
| Upholstery nozzle                          | •                      |
| Crevice nozzle                             | •                      |
| Dusting brush with synthetic bristles      | •                      |
| Floorhead SBD 365-3                        | •                      |
| Turbobrush STB 305-3 (DOUBLE; LOCKED)      | •                      |
| Vacuum cleaner bag voucher                 | •                      |