

Guard M1  
Cylinder vacuum cleaner  
powerful suction | AirClean filter | integrated accessories | height-adjustable



- ComfortFit: developed for the highest hygiene standards
- Durable: tested for the equivalent of 20 years\*\* use
- 99.99 %\*\* dust retention capacity with AirClean filter
- Mode button: simply adjust the suction power to your needs
- Practical compartment with three-piece accessory set for added convenience
- Stylish design concept and illuminated LED icons
- Range extra, fewer socket changes: 11 m operating radius

\* Compared to SBB300-3 Parquet Twister parquet brush

\*\* In accordance with IEC 62885-2:2021

EAN: 4002516832942 / Material number: 12652510 /  
Old Material Number: 41VWF030GB

<b>Construction type and design</b>	
Colour	Terra red
<b>Construction type</b>	
Cylinder vacuum cleaner with dustbag	•
Product brand	Guard M1
Motor type	PowerLine
Product type	SVWFO
<b>User convenience</b>	
Electronic suction power control	MODE button
Ergonomic carrying handle	•
Ergonomic handle	•
Suction tube	HES mid-size telescopic tube
Park system for pauses during vacuuming	•
Park system	on both sides
Comfort cable rewind	•
<b>Mobility</b>	
Castors with steel shaft, gentle on floors	•
3-point manoeuvrability	•
Operating radius in m	11
Cable length in m	7.5
<b>Filter system</b>	
Vacuum cleaner bag type	HyClean Pure CO
Dustbag volume in l	4.5
Dustbag change indicator	•
Exhaust filter	AirClean filter
Motor protection filter	•
<b>Safety</b>	
Locking system for floorhead and handle	•
Automatic dustbag positioning	•
All-round furniture protection strip	•
<b>Technical data</b>	
Power	890
Volume in dB	77
Width in mm	496
Height in mm	227
Weight in kg	6.7
<b>Accessories included</b>	
Three-piece accessory set	Integrated
Upholstery nozzle	•
Crevice nozzle	•
Dusting brush with synthetic bristles	•
Floorhead	•