

Guard M1
Cylinder vacuum cleaner
powerful suction | AirClean filter | integrated accessories | height-adjustable



- ComfortFit: developed for the highest hygiene standards
- Durable: tested for the equivalent of 20 years** use
- 99.99 %** dust retention capacity with AirClean filter
- Mode button: simply adjust the suction power to your needs
- Practical compartment with three-piece accessory set for added convenience
- Stylish design concept and illuminated LED icons
- Range extra, fewer socket changes: 11 m operating radius

* During the development phase of the advertised Guard vacuum cleaner series, Miele tested models and core components at various power levels for 750 hours of use (= 37.5 hours p.a.). More:

miele.com/20years

** In accordance with IEC 62885-2:2021

EAN: 4002516832942 / Material number: 12652510 /
Old Material Number: 41VWF030GB

| | |
|--|------------------------------|
| General information | |
| Model identifier | Guard M1 |
| Manufacturer material number | 12652510 |
| EAN number | 4002516832942 |
| Construction type | |
| Cylinder vacuum cleaner with dustbag | • |
| Product brand | Guard M1 |
| Motor type | PowerLine |
| Product type | SVWFO |
| Construction type and design | |
| Colour | Terra red |
| User convenience | |
| Electronic suction power control | MODE button |
| Ergonomic carrying handle | • |
| Ergonomic handle | • |
| Suction tube | HES mid-size telescopic tube |
| Park system for pauses during vacuuming | • |
| Park system for storing | on both sides |
| Comfort cable rewind | • |
| Mobility | |
| Castors with steel shaft, gentle on floors | • |
| 3-point manoeuvrability | • |
| Operating radius in m | 11 |
| Cable length in m | 7.5 |
| Filter system | |
| Vacuum cleaner bag type | HyClean Pure CO |
| Dustbag volume in l | 4.5 |
| Dustbag change indicator | • |
| Exhaust filter | AirClean filter |
| Motor protection filter | • |
| Safety | |
| Locking system for floorhead and handle | • |
| Automatic dustbag positioning | • |
| Furniture protection | • |
| Technical data | |
| Power | 890 |
| Volume in dB | 77 |
| Width in mm | 496 |
| Height in mm | 227 |
| Depth in mm | 284 |
| Weight in kg | 6.7 |
| Accessories included | |
| Accessory storage | Integrated |
| Upholstery nozzle | • |
| Crevice nozzle | • |
| Dusting brush | • |
| Floorhead SBD 365-3 | • |

Guard M1

Cylinder vacuum cleaner

powerful suction | AirClean filter | integrated accessories | height-adjustable

EAN: 4002516832942 / Material number: 12652510 /

Old Material Number: 41VWF030GB

HyClean Pure CO vacuum cleaner bag

Voucher for HyClean Pure CO

•

•

