

Guard S1
Cylinder vacuum cleaner
high suction | light, compact | in modern trend colours | easy to adjust height



- ComfortFit: developed for the highest hygiene standards
- Ergonomic handle: for comfortable handling
- 99.99 %* dust retention capacity with AirClean filter
- Durable: tested for the equivalent of 20 years*** use
- 10 m operating radius: effortless cleaning of larger areas
- Easy handling thanks to lightweight and compact design
- Height-adjustable EasySlide telescopic tube

* In accordance with IEC 62885-2:2021

** During the development phase of the advertised Guard vacuum cleaner series, Miele tested models and core components at various power levels for 750 hours of use (= 37.5 hours p.a.). More:

miele.com/20years

EAN: 4002516833741 / Material number: 12652890 /
Old Material Number: 41TAF032GB

General information	
Model identifier	Guard S1
Manufacturer material number	12652890
EAN number	4002516833741
Construction type	
Cylinder vacuum cleaner with dustbag	•
Product brand	Guard S1
Motor type	PowerLine
Product type	STAFO
Appliance type and design	
Appliance colour	Nordic blue
User convenience	
Electronic suction power control	4-stage rotary selector
Suction tube	EasySlide telescopic tube
Park system for pauses during vacuuming	•
Mobility	
Castors with steel shaft, gentle on floors	•
3-point manoeuvrability	•
Operating radius in m	10
Cable length in m	7.5
Filter system	
Vacuum cleaner bag type	HyClean Pure TU
Dustbag volume in l	3.5
Dustbag change indicator	•
Exhaust filter	AirClean filter
Motor protection filter	•
Safety	
Locking system for floorhead and handle	•
Automatic dustbag positioning	•
Furniture protection	•
Technical data	
Power	890
Volume in dB	76.0
Width in mm	425
Height in mm	230
Depth in mm	245
Weight in kg	6.0
Accessories included	
Upholstery nozzle	•
Crevice nozzle	•
Dusting brush	•
Floorhead SBD 355-3	•
HyClean Pure TU vacuum cleaner bag	•
Voucher for HyClean Pure TU	•