


U 504

Basket

For the optimum loading of cutlery.



- With 9 plastic cutlery holders
- Compartment size 106 x 106 mm

EAN: 4002516854593 / Material number: 12655860 /
Old Material Number: 67650402D

Product category

Fresh-water dishwashers

•

Tank dishwashers

•

Throughfeed tank dishwashers

•

Machine affiliation

G 8066

•

PG 8041

•

PG 8043

•

PG 8045

•

PG 8055

•

PG 8055 U

•

PG 8056

•

PG 8056 U

•

PG 8057 TD

•

PG 8057 TD U

•

PG 8058

•

PG 8058 U

•

PG 8059

•

PG 8059 U

•

PG 8060

•

PG 8061

•

PG 8061 U

•

PG 8165

•

PG 8166

•

PG 8169

•

PG 8172

•

PG 8172 Eco

•

PTD 901

•

PG 8096

•

PG 8096 U

•

PG 8099

•

PG 8099 U

•

PFD 400

•

PFD 400 U

•

PFD 401

•

PFD 401 U

•

PFD 402

•

PFD 402 U

•

PFD 404

•

PFD 404 U

•

PFD 405

•

PFD 405 U

•

PFD 407

•

PFD 407 U

•

PFD 408 U

•

PTD 702

•

PTD 703

•

PTD 704

•

U 504

Basket

For the optimum loading of cutlery.



EAN: 4002516854593 / Material number: 12655860 /
Old Material Number: 67650402D

Product category

Basket for cutlery

•

Product properties

Material	Wire, rilsan-coated, Plastic
Colour	Grey, White
Number of compartments	9
Compartment width in mm	106
Compartment length in mm	106

Dimensions and weight

External dimensions, gross height in mm	180
External dimensions, gross width in mm	520
External dimensions, gross depth in mm	610
Net weight in kg	2.27
Gross weight in kg	3.29