

U 634-S

Basket

with tilters for the improved washing and drying of 16 glasses.

EAN: 4002516855217 / Material number: 12663010 / Old Material Number: 67663402D



- 16 compartments, compartment size 113 mm x 113 mm mm
- Fitted with a low top frame
- Height 104/173/233, width 500 (incl. GN sides 530), depth 500 mm
- Includes extender set comprising short and long extenders

Product category	
Fresh-water dishwashers	•
Tank dishwashers	•
Throughfeed tank dishwashers	•
Machine affiliation	
G 8066	•
PG 8067	•
PG 8041	•
PG 8043	•
PG 8045	•
PG 8055	•
PG 8055 U	•
PG 8056	•
PG 8056 U	•
PG 8057 TD	•
PG 8057 TD U	•
PG 8058	•
PG 8058 U	•
PG 8059	•
PG 8059 U	•
PG 8060	•
PG 8061	•
PG 8061 U	•
PG 8165	•
PG 8166	•
PG 8169	•
PG 8172	•
PG 8172 Eco	•
PTD 901	•
PG 8096	•
PG 8096 U	•
PG 8099	•
PG 8099 U	•
PFD 400	•
PFD 400 U	•
PFD 401	•
PFD 401 U	•
PFD 402	•
PFD 402 U	•
PFD 404	•
PFD 404 U	•
PFD 405	•
PFD 405 U	•
PFD 407	•
PFD 407 U	•
PFD 408 U	•
PTD 702	•
PTD 703	•

U 634-S

Basket

with tilters for the improved washing and drying of 16 glasses.



EAN: 4002516855217 / Material number: 12663010 / Old Material Number: 67663402D

PTD 704

•

Product category

Basket for glassware

•

Product properties

Material	Plastic
Colour	Blue, Yellow
Number of compartments	16
Compartment width in mm	113
Compartment length in mm	113

Capacity

Glasses 110–140 mm high [number]	16
Glasses 175–200 mm high [number]	16

Dimensions and weight

External dimensions, gross height in mm	220
External dimensions, gross width in mm	550
External dimensions, gross depth in mm	550
Net weight in kg	2.72
Gross weight in kg	3.66