


E 167  
Cutlery basket  
For holding cutlery.



- Suitable for plate basket U 506
- 7 Compartments

EAN: 4002516854869 / Material number: 12663060 / Old Material Number: 67616702D

Product category	
Fresh-water dishwashers	•
Tank dishwashers	•
Throughfeed tank dishwashers	•
Machine affiliation	
G 8066	•
PG 8041	•
PG 8043	•
PG 8045	•
PG 8055	•
PG 8055 U	•
PG 8056	•
PG 8056 U	•
PG 8057 TD	•
PG 8057 TD U	•
PG 8058	•
PG 8058 U	•
PG 8059	•
PG 8059 U	•
PG 8060	•
PG 8061	•
PG 8061 U	•
PG 8165	•
PG 8166	•
PG 8169	•
PG 8172	•
PG 8172 Eco	•
PTD 901	•
PG 8096	•
PG 8096 U	•
PG 8099	•
PG 8099 U	•
PFD 400	•
PFD 400 U	•
PFD 401	•
PFD 401 U	•
PFD 402	•
PFD 402 U	•
PFD 404	•
PFD 404 U	•
PFD 405	•
PFD 405 U	•
PFD 407	•
PFD 407 U	•
PFD 408 U	•
PTD 702	•
PTD 703	•
PTD 704	•

E 167  
Cutlery basket  
For holding cutlery.



EAN: 4002516854869 / Material number: 12663060 / Old Material Number: 67616702D

<b>Product category</b>	
Insert for cutlery	•
<b>Product properties</b>	
Material	Plastic
Colour	White
Number of compartments	7
<b>Dimensions and weight</b>	
External dimensions, gross height in mm	95
External dimensions, gross width in mm	170
External dimensions, gross depth in mm	550
Net weight in kg	0.19
Gross weight in kg	0.34